

SET-UP & POSITIONER Rx

Lab use only
POS
SU
MT LOW
DUP
BA
IMP
CAR
HA
SER
PC
PL
SC

Dr. _____
 Address _____
 _____ Zip _____
 Patient _____
 Date Shipped to PAR _____
 Date Required _____

- Standard Set-up and Positioner
- Articulator Set-up and Positioner
- Positioner Only
- Set-up Only
- Pre-treatment Diagnostic Set-up
- Other _____
- Custom Athletic Mouthguard for:

_____ (name of sport)

- One Piece upper and lower
- Single arch vacuum form

SET-UP INSTRUCTIONS:

- Reset teeth on enclosed models
- Duplicate our models for setup
- Return Controls
- Do not carve brackets and bands
- Allow for lingual 3-3 arch on set-up
- Reset only teeth marked:



ARCH WIDTH:

- Upper: Maintain Constrict Expand
 Lower: Maintain Constrict Expand

SPACE CLOSURE:

- Close all completely
- Close as feasible
- Leave space between _____

ANTERIOR ROOT INCLINATION:

- Upper: Maintain Lingual _____ ° Labial _____ °
 Lower: Maintain Lingual _____ ° Labial _____ °

OCCLUSAL PLANE:


- Flat Curve of Spee As is

VERTICAL OVERBITE:

- As is Reduce _____ mm Increase _____ mm

HORIZONTAL OVERJET:

- As is Reduce Increase



SPECIAL INSTRUCTIONS

POSITIONER INSTRUCTIONS:

POSITIONER OR MOUTHGUARD MATERIAL:

- FLEXICLEAR:** (Unflavored) Soft Medium Firm
 Available with: Mint
 Colors with Flavors: Available in medium Flexiclear only.
Call for Colors

- Impak - clear (heat to soften)
- Silicone - translucent, elastic

TRIM: Height

- High (5mm above gingival)
- Standard (1-3 mm above gingival)
- Short (to gingival)
- Mini-trim (3/4 of crowns)

- Thickness: Standard Thick Thin

END POSITIONER DISTAL TO:

- 1st Molars 2nd Molars Other _____
- (NOTE: Positioner should cover all teeth to prevent super-eruption.)

ARTICULATION:

- Average bite opening
- Hinge Axis Tracing Headfilm
- Core Technique
- Adjustable Articulator - indicate below:
- Denar Sam Whip Mix Panadent Hanau

OPTIONS:

- Air Holes
- PAR Clasps (for Placement And Retention)

Mark position of clasps

between teeth on diagram: R 7 6 5 | 5 6 7 L

- Ball End Loop End

- PAR Socket Liners on teeth marked: R 2 1 | 1 2 L

- Posterior Serrations

PLEASE SEND ADDITIONAL SUPPLIES:

- Boxes Rx Forms Other